



## **User Manual**

[Revision 1.0]

---

**READ THIS MANUAL CAREFULLY BEFORE USE – FAILURE TO DO SO MAY RESULT IN INJURY, PROPERTY DAMAGE AND MAY VOID WARRANTY. • KEEP THIS MANUAL FOR FUTURE REFERENCE. • Products covered by this manual may vary in appearance, assembly, inclusions, specifications, description, and packaging.**

---

# Safety

Safety messages are designed to alert you to possible dangers or hazards that could cause death, injury or equipment or property damage if not understood or followed. Safety messages have the following symbols:



You **WILL** be **KILLED** or **SERIOUSLY INJURED** if you do not follow instructions.



You **CAN** be **KILLED** or **SERIOUSLY INJURED** if you do not follow instructions.



You **CAN** be **INJURED** if you do not follow instructions or equipment damage may occur.

It is vital that you read and understand this user manual before using the product, including safety warnings, and any assembly and operating instructions. Keep the manual for future reference.

Safety precautions and recommendations detailed here must be fully understood and followed to reduce the risk of injury, fire, explosion, electrical hazard, and/or property damage.

Safety information presented here is generic in nature – some advice may not be applicable to every product. The term "equipment" refers to the product, be it electrical mains powered, battery powered or combustion engine powered.

- **Before Use** - If you are not familiar with the safe operation/handling of the equipment or are in any way unsure of any aspect of suitability or correct use for your application, you should complete training conducted by a person or organization qualified in safe use and operation of this equipment, including fuel/electrical handling and safety.
- Do NOT operate the equipment in flammable or explosive environments, such as in the presence of flammable liquids, gases, or dust. The equipment may create sparks or heat that may ignite flammable substances.
- Keep clear of moving parts.
- Equipment may be a potential source of electric shock or injury if misused.
- Do NOT operate the equipment if it is damaged, malfunctioning or is in an excessively worn state.
- Do NOT allow others to use the equipment unless they have read this manual and are adequately trained.
- Keep packaging away from children - risk of suffocation! Operators must use the equipment correctly. When using the equipment, consider conditions and pay due care to persons and property.

### General Work Area Safety

- Work areas should be clean and well lit.
- Do not operate the equipment if bystanders, animals etc are within operating range of the equipment or the general work area.
- If devices are provided for connecting dust extraction / collection facilities, ensure these are connected and used properly. Dust collection can reduce dust-related hazards.

### General Personal Safety

- Wear appropriate protective equipment when operating, servicing, or when in the operating area of the equipment to help protect from eye and ear injury, poisoning, burns, cutting and crush injuries. Protective equipment such as safety goggles, respirators, non-slip safety footwear, hard hat, hearing protection etc should be used for appropriate equipment / conditions. Other people nearby should also wear appropriate personal protective equipment. Do not wear loose clothing or jewellery, which can be caught in moving parts. Keep hair and clothing away from the equipment.
- Stay alert and use common sense when operating the equipment. Do not over-reach. Always maintain secure footing and balance.
- Do not use the equipment if tired or under the influence of drugs, alcohol, or medication.
- This equipment is not intended for use by persons with reduced physical, sensory, or mental capabilities.

### General Fuel Safety

- Petrol/fuel/gasoline is extremely flammable – keep clear of naked flames or other ignition sources.
- Do not spill fuel. If you spill fuel, wipe it off the equipment immediately – if fuel gets on your clothing, change clothing.
- Do NOT smoke near fuel or when refuelling.
- Always shut off the engine before refuelling.
- Do NOT refuel a hot engine.
- Open the fuel cap carefully to allow any pressure build-up in the tank to release slowly.
- Always refuel in well-ventilated areas.
- Always check for fuel leakage. If fuel leakage is found, do not start, or run the engine until all leaks are fixed.

### General Carbon-Monoxide Safety

- Using a combustion engine indoors **CAN KILL IN MINUTES**. Engine exhaust contains carbon-monoxide – a poison you cannot smell or see.
- Use combustion engines **OUTSIDE** only, and far away from windows, doors, and vents.

### General Equipment Use and Care

- The equipment is designed for domestic use only.
- Handle the equipment safely and carefully.
- Before use, inspect the equipment for misalignment or binding of moving parts, loose components, damage, or any other condition that may affect its operation. If damaged, have the equipment repaired by an authorised service centre or technician before use.
- Prevent unintentional starting of the equipment - ensure equipment and power switches are in the OFF position before connecting or moving equipment. Do not carry equipment with hands or fingers touching any controls. Remove any tools or other items that are not a part of the equipment from it before starting or switching on.
- Do not force the equipment. Use the correct equipment for your application. Equipment will perform better and be safer when used within its design and usage parameters.
- Use the equipment and accessories etc. in accordance with these instructions, considering working conditions and the work to be performed. Using the equipment for operations different from those intended could result in hazardous situations.
- Always keep equipment components (engines, hoses, handles, controls, frames, housings, guards etc) and accessories (cutting tools, nozzles, bits etc) properly maintained. Keep the equipment clean and, where applicable, properly lubricated.
- Store the equipment out of reach of children or untrained persons. To avoid burns or fire hazards, let the equipment cool completely before transporting or storing. Never place or store the equipment near flammable materials, combustible gases, or liquids etc.
- The equipment is not weather-proof, and should not be stored in direct sunlight, at high ambient temperatures or locations that are damp or humid.
- Do not clean equipment with solvents, flammable liquids, or harsh abrasives.
- For specific equipment safety use and care, see Equipment Safety.

**General Electrical Safety**

- Inspect electrical equipment, extension cords, power bars, and electrical fittings for damage or wear before each use. Repair or replace damaged equipment immediately.
- Ensure all power sources conform to equipment voltage requirements and are disconnected before connecting or disconnecting equipment.
- When wiring electrically powered equipment, follow all electrical and safety codes.
- Always use a Residual Current Device (RCD) (GFCI) / Earth Leakage Circuit Breaker / Safety Switch when operating this item (electrical items only).
- High voltage / high current power lines may be present. Use extreme caution to avoid contact or interference with power lines. Electrical shock can be fatal.

**General Electrical Safety**

- Electrically grounded equipment must have an approved cord and plug and be connected to a grounded electrical outlet.
- Do NOT bypass the ON/OFF switch and operate equipment by connecting and disconnecting the electrical cord.
- Do NOT use equipment that has exposed wiring, damaged switches, covers or guards.
- Do NOT use electrical equipment in wet conditions or in damp locations.
- Do NOT use electrical cords to lift, move or carry equipment.
- Do NOT coil or knot electrical cords, and ensure electrical cords are not trip hazards.

**General Service Information**

- The equipment must be serviced or repaired at authorised service centres by qualified personnel only.
- Replacement parts must be original equipment manufacturer (OEM) to ensure equipment safety is maintained.
- Do NOT attempt any maintenance or repair work not described in this manual.
- After use, the equipment and components may still be hot – allow the equipment to cool and disconnect spark plugs and/or electrical power sources and/or batteries from it before adjusting, changing accessories, or performing repair or maintenance.
- Do NOT adjust while the equipment is running.
- Perform service-related activities in suitable conditions, such as a workshop.
- Replace worn, damaged or missing warning/safety labels immediately.

**Child Choking Hazard!**

**• WARNING: CHOKING HAZARD - CHILDREN UNDER 8 YRS. CAN CHOKE OR SUFFOCATE ON SMALL PARTS. ADULT SUPERVISION REQUIRED.**

- Carefully inspect anything a baby or child could gain access to, including this product & packaging.
- Dangers include anything of a size that could become stuck or lodged in a baby or child's airway.
- Dangers include things like raw, hard fruit and vegetable pieces, large pieces of meat, bones or sausage skins, popcorn, nuts, hard lollies and corn chips, small magnets and batteries, coins, beads, marbles and small uninflated balloons, broken toys and smaller toys, buttons, keys, and coins. But there are also many more.
- Pay particular attention to accessories, nuts, bolts, screws, washers, caps, covers, loose parts, parts that could be removed or break off, pieces of packaging, staples, tape etc.
- In the case of any doubt, secure the item of concern in an appropriate manner, or completely remove the possibility of access by a baby or child.

**Important Transfer Chair Safety Information**

- When attempting to stand you must ensure the footrests (when equipped) have been flipped up prior to attempting to exit the chair.
- Assistance may be required to exit the chair. A qualified carer is recommended when traveling up or down an incline.
- Do not stand directly on the wheels, do not place any weight directly on the wheels.
- When sitting on the seat ensure the brakes are engaged and the user is balanced towards the back of the seat. Do not lean forward.
- The rollator is not designed for transiting a patient through long distances or for use on rough, uneven surfaces.

**Additional Safety Reminders**





















- This product is intended for the disabled, the sick, and the elderly with mobility difficulties, for transfer or urination purposes.
- Be sure to read these "Safe Use Precautions" carefully before use to ensure correct operation.
- Failure to follow the instructions in "Safe Use Precautions" may result in product damage, injury, or even life-threatening situations.
- The operator should understand the use of the manual before use.
- Do not allow young children or individuals who cannot understand the operation to use the transfer chair.
- The "★" symbol refers to non-standard configurations. Each of these precautions is important for safety, so be sure to follow them.
- After reading the instruction manual, keep it in a place where it can be easily accessed at all times.
- Pre-use Check: Ensure all components are complete and securely installed for safe use.
- Height Adjustment: Adjust the chair to the appropriate height before use, ensuring both sides of the seat plate remain at the same height.
- Entering and Exiting: Do not stand on the footboard when the user sits in or leaves the transfer chair. Step on the ground to prevent the transfer chair from tipping over.
- Safety Latch and Seat Belt: Ensure the safety latch and seat belt are locked when using the transfer chair.
- Wheel Locking: Ensure the wheels are locked when the user is entering or leaving the transfer chair.
- Reaching for Items: When the user sits in the transfer chair and wants to take items outside the chair (in front, on both sides, or behind), avoid leaning too far to prevent tipping.
- Picking Up Objects: Do not bend below your knees to pick up objects from the floor.
- Seat Height Adjustment: Lock the front wheel when adjusting the seat height and move both handles uniformly at the same speed.
- Avoid Inertia: Do not use inertia to rush the transfer chair up or down steps to prevent accidents or damage. Do not move the user together with the transfer chair when going up or down stairs to avoid danger.













## Manual Transfer Chair

- **Movement:** The transfer chair is suitable for moving on flat ground. Avoid potholes and obstacles to prevent tipping.
- **Downhill Slope:** When the downhill slope is greater than 5 degrees, turn the transfer chair backward and carefully proceed. Facing forward on a downhill slope may cause the user to lean forward and fall.
- **Staying on Sloping Ground:** Lock all four wheels when staying on sloping ground.
- **Brake Usage:** The brakes on the four wheels of the transfer chair are only for docking and not for deceleration or braking during travel.

# Safety Symbols

The product may have safety warning labels attached to it, explained below. Understand the symbols on your product and their meanings. If any stickers become unreadable, unattached etc., replace them.

 <p><b>Flammable Material Hazard</b> Flammable liquids, gases, or substances etc may present. Avoid ignition sources and open flames. Danger of fire.</p>	 <p><b>Read User Manual</b> Read and fully understand product safety warnings, operation, procedures etc before using the product.</p>	 <p><b>Use Hand Protection</b> Wear appropriate hand protection and take due care as the product or use of the product may present hand hazards.</p>	 <p><b>Carbon-Monoxide Hazard</b> Do not use the product in confined areas or without adequate ventilation. Carbon-monoxide poisoning can be fatal.</p>
 <p><b>Electrocution / Electrical Shock Hazard</b> High voltage or high current electricity may be present or required by the product. Take due care when handling electrical products, cables, plugs and leads. Electrical shock can be fatal.</p>	 <p><b>Toxic Fumes / Dust Hazard</b> Using the product or by-products from use may produce fumes, smoke or particles that could be harmful if inhaled. Wear appropriate breathing protection and have adequate ventilation.</p>	 <p><b>Explosive Material Hazard</b> Combustible liquids, gases, or substances etc may be present. Avoid ignition sources and open flames. Danger of explosion.</p>	 <p><b>Cutting / Amputation Hazard</b> The product may have blades, edges or mechanical devices that can cause severe cut injury to fingers, limbs etc. Take due care when handling and using the product.</p>
 <p><b>Crush Hazard</b> The product may have blades, edges or mechanical devices that can cause severe crush injury to fingers, limbs etc. Take due care when handling and using the product.</p>	 <p><b>Single Operator Only</b> The product must be operated by a single person only. More than one person operating the product may introduce additional hazards.</p>	 <p><b>Use Face Protection</b> Wear appropriate full-face protection and take due care as the product or use of the product may present face and eye hazards.</p>	 <p><b>Use Foot Protection</b> Wear appropriate foot protection and take due care as the product or use of the product may present foot hazards.</p>
 <p><b>Use Eye / Ear / Head Protection</b> Wear appropriate eye and / or ear and / or head protection and take due care as the product or use of the product may present eye, hearing, and head hazards.</p>	 <p><b>Running Hazard</b> Do not run on or near the product as doing so may present a fall hazard.</p>	 <p><b>Diving Hazard</b> Do not dive into the product as doing so may present a neck / head injury hazard.</p>	 <p><b>Adult Supervision Required</b> Always supervise children and other users of a product to prevent drowning or injury.</p>
 <p><b>Skin Penetration / Puncture Hazard</b> The product may produce pressure, emit liquids or objects that can cause severe injury to fingers, limbs, blood etc. Take due care when handling and using the product.</p>	 <p><b>Hot Surface Hazard</b> Be aware that the product may produce high temperatures and hot surfaces that can cause burn injuries.</p>	 <p><b>Flying Debris Hazard</b> Be aware that the product or use of the product may present hazards produced by flying debris. Wear appropriate clothing and protective devices.</p>	 <p><b>Moving Parts Hazard</b> Be aware that the product contains or uses mechanical devices that move or rotate. Always wait for moving parts to stop fully before handling the product, adjusting, maintenance etc.</p>

 <p><b>Carbon-Monoxide Hazard</b> Do not use the product in confined areas or without adequate ventilation. Carbon-monoxide poisoning can be fatal.</p>	 <p><b>Pull Hazard</b> Be aware that the product contains or uses mechanical devices that can pull in objects and can cause severe injury to fingers, limbs etc. Take due care when handling and using the product.</p>	 <p><b>Slope / Fall Injury Hazard</b> Be aware that using the product on sloping surfaces or in slippery conditions may present additional dangers from falls and contact with blades, moving parts, hot surfaces etc.</p>	 <p><b>"Slam Dunk" Warning</b> Do NOT attempt "slam dunk" manoeuvres as this may result in severe injury due to falling, product breakage or collapse etc.</p>
 <p><b>Electrocution / Electrical Shock Hazard - Outdoor</b> High voltage or high current electricity may be present or required by the product. Do NOT use in rain, damp, or wet conditions. Electrical shock can be fatal.</p>	 <p><b>Electrocution / Electrical Shock Hazard - Disconnect</b> High voltage or high current electricity may be present or required by the product. Always disconnect the product from the electrical supply before handling the product, adjusting, maintenance etc.</p>	 <p><b>Power Line Electrocution Hazard</b> High voltage / high current power lines may be present. Use extreme caution to avoid contact or interference with power lines. Electrical shock can be fatal.</p>	 <p><b>"Kick-Back" Hazard</b> High level of "kick-back" hazard that can cause the machine to suddenly rotate towards operator. Kick-back injury can be fatal.</p>
 <p><b>Winch Operator Position Hazard</b> Do NOT stand between winch and load. Do NOT use winch to move people.</p>	 <p><b>Winch Lift Hazard</b> Do NOT LIFT load vertically. Use machine to PULL only.</p>	 <p><b>Cable Hazard</b> Ensure that load bearing cable is not kinked or knotted.</p>	 <p><b>Winch Cable Hazard</b> Ensure that there is a minimum number of cable coils on winching mechanism.</p>
 <p><b>Winch Hook Hazard</b> Carry hook to load – do NOT throw or run.</p>	 <p><b>Flash / Blinding Hazard</b> Wear appropriate eye protection for welding. Direct exposure to weld arcs may cause permanent eye injury.</p>	 <p><b>Laser Hazard</b> Laser may be in use – do NOT look directly at laser or allow others to.</p>	 <p><b>Choking Hazard</b> Children under 8 yrs. Can choke or suffocate on small parts. Adult supervision required.</p>

# Table of Contents

<b>Safety</b> .....	<b>2</b>
Safety Symbols .....	5
<b>Parts Identification</b> .....	<b>8</b>
<b>Assembly</b> .....	<b>9</b>
<b>Operation</b> .....	<b>11</b>
<b>Maintenance and Storage</b> .....	<b>15</b>
<b>Specifications</b> .....	<b>16</b>

# Parts Identification




No.	Name
1	Front wheels
2	Footboard
3	Rear wheel tube
4	Rear wheels
5	Base
6	Commode / Toilet bucket
7	Seat cushion
8	Left seat frame
9	Right seat frame

# Assembly

1. Cut the cable tie (FIG. 1) fixing the pinion.
2. Turn the pinion downward and insert the D-shaped pin (FIG. 2) through the base from outside to inside the pin hole.
3. Insert the pin into the pin to prevent it from falling off (FIG. 3).
4. Install the pinion in the same way on the other side.



(FIG. 1)



(FIG. 2)

5. Hold the seat frame by hand, align the rotating sleeve of the seat frame with the rotating support sleeve of the telescopic frame, and press the cap end of the quick-disconnecting shaft downward with your thumb and insert it (Figure 4). Ensure that the cap end of the quick-disconnecting shaft is pressed after insertion. The seat frame should not be removed when lifted upward if properly installed.




(FIG. 3)



(FIG. 4)

6. If you need to use a toilet bucket, first install the toilet bucket wire rack into the mounting holes of the two columns. Then, load the toilet bucket from the rear, as shown in the following figure:



★ If the infusion stand is optional, it can be installed on the outside of the left or right seat.

7. After all the accessories are installed, turn the handle in the opposite direction, lift the handle up, and you can push the transfer chair normally.



# Operation

**IMPORTANT!** When attempting to stand you must ensure the footrests (when equipped) have been flipped up prior to attempting to exit the chair. Assistance may be required to exit the chair. A qualified carer is recommended when traveling up or down an incline. Do not stand directly on the wheels, do not place any weight directly on the wheels. When sitting on the seat ensure the brakes are engaged and the user is balanced towards the back of the seat. Do not lean forward. The rollator is not designed for transiting a patient through long distances or for use on rough, uneven surfaces.

The following steps explain the operation for transferring a user from the seat to the bed. The method for transferring between other furniture and facilities is the same.

1. Ensure the user is sitting upright. Adjust the height of the transfer chair frame so that the bottom of the frame is slightly higher than the seat plane.
2. Open both sides of the seat and move the chair to the front of the user. Adjust the position of the transfer chair so that the user's centre of gravity is in the centre of the seat. Adjust the height of the seat again so that the bottom of the seat slightly touches the seat plane.
3. Lock the 5-inch front wheel brake, lift the user's feet, and place them on the foot pedal.



**Manual Transfer Chair**

4. Let the user tilt their body slightly to the left and right in turn, while the operator rotates the corresponding side of the seat frame. This allows the user to sit on the left and right seat frames of the transfer chair in turn.
5. When the user is fully seated on the seat frame and the left and right seat frames are completely closed, the safety lock will automatically engage. The operator should then lock the seat belt.
6. The user should adjust to a comfortable sitting position and grasp the handrails on both sides with both hands.
7. The operator checks the distance between the bottom of the transfer chair frame and the seat. Raise the seat frame by 2-3 cm, release the 5-inch front wheel brake, and slowly move the transfer chair away from the seat.



8. The operator pushes the transfer chair from behind to the edge of the bed.
9. Check the height of the transfer chair at the edge of the bed to ensure that the bottom of the seat is slightly higher than the bed. If the height of the seat needs to be adjusted, the front wheel must be locked before adjustment.
10. The operator pushes the transfer chair completely to the predetermined position by the bed, locks the 5-inch front wheel brake, and adjusts the height of the seat again so that the bottom of the seat slightly touches the bed.
11. Open the seat belt and safety lock on the back of the transfer chair, slowly open the left and right seat frames, and allow the user's buttocks to sit on the bed.



12. The operator releases the 5-inch front wheel lock, lifts the user's feet, and moves the transfer chair away from the bed.

★ If the infusion rod is selected, loosen the knob to adjust the height of the infusion rod. When the infusion rod is installed in the transfer chair, it needs to be removed if the height of the seat needs to be adjusted.

★ If the safety belt is optional, adjust the length to the appropriate position before use.



# Maintenance and Storage

## Maintenance:

- Do not scratch or scrape the surface of the product with hard objects to prevent material damage.
- For daily cleaning and maintenance, use a neutral synthetic detergent. Wipe with a wet cloth and then dry with a dry cloth.
- Avoid contact with strong acids and alkalis to prevent corrosion.
- Regularly check for loose nuts and screws, as well as damaged or broken parts. Replace broken parts if necessary.

## Storage:

Packaged products should be stored in a temperature range of 0°C to 40°C, with relative humidity not exceeding 80%. Store in a well-ventilated, clean room free of corrosive gases. During transportation, avoid exposure to moisture, light, and heat sources.

# Specifications

<b>Type</b>	Manual
<b>Colour</b>	White, Grey
<b>Seat Width</b>	46cm
<b>Seat Depth</b>	40cm
<b>Max User Load</b>	120kg
<b>Frame material</b>	2mm double thickened steel
<b>Recommended User Height</b>	135cm - 198cm
<b>Rear wheels</b>	3" 360° castor wheels
<b>Front wheels</b>	5" castor wheels w/ brakes
<b>Seat Height Range</b>	39 - 59cm
<b>Height of base</b>	12cm
<b>IP rating</b>	IPX5 (protected against low pressure water stream)

*Note: Requires minimal assembly*



**Some experts believe that the incorrect or prolonged use of almost any product may cause serious injury or death. To help reduce your risk of serious injury or death, refer to the information below. For more information, see [www.datastreamserver.com/safety](http://www.datastreamserver.com/safety)**

- Consult all documentation, packaging, and product labelling before use. Note that some products feature documentation available online. It is recommended to print and retain the documentation.
- Before each use, check the product for loose/broken/damaged/missing parts, wear, or leaks (if applicable). Never use a product with loose/broken/damaged/missing parts, wear, or leaks.
- Products must be inspected and serviced (if applicable) by a qualified technician every 6 months. This is based on average residential use by persons of average size and strength, and on a property of average metropolitan size. Use beyond these recommendations may require more frequent inspections/servicing.
- Ensure that all users of the product have completed a suitable industry recognised training course before being allowed access to the product.
- The product has been supplied by a general merchandise retailer that may not be familiar with your specific application or description of application. Be sure to attain third-party approval from a qualified specialist for your application before use, regardless of any assurances from the retailer or its representatives.
- This product is not intended for use where fail-safe operation is required. As with any product (for example, automobile, computer, toaster), there is the possibility of technical issues that may require the repair or replacement of parts, or the product itself. If the possibility of such failure and the associated time it may take to rectify could in any way inconvenience the user, business, or employee, or financially affect the user, business, or employee, then the product is not suitable for your requirements. This product is not intended for use where incorrect operation or a failure of any kind, including but not limited to, a condition requiring product return, replacement, parts replacement, or service by a technician may cause financial loss, loss of employee time or an inconvenience requiring compensation.
- If this product has been purchased in error when considering the information presented here, contact the retailer directly for details of their returns policy, if required.

