

Premium Pediatric Wheelchair

EQUIPMED®



User Manual

[Revision 2.0]

READ THIS MANUAL CAREFULLY BEFORE USE – FAILURE TO DO SO MAY RESULT IN INJURY, PROPERTY DAMAGE AND MAY VOID WARRANTY. • KEEP THIS MANUAL FOR FUTURE REFERENCE. • Products covered by this manual may vary in appearance, assembly, inclusions, specifications, description, and packaging.

Safety

Safety messages are designed to alert you to possible dangers or hazards that could cause death, injury or equipment or property damage if not understood or followed. Safety messages have the following symbols:



You **WILL** be **KILLED** or **SERIOUSLY INJURED** if you do not follow instructions.



You **CAN** be **KILLED** or **SERIOUSLY INJURED** if you do not follow instructions.



You **CAN** be **INJURED** if you do not follow instructions or equipment damage may occur.

It is vital that you read and understand this user manual before using the product, including safety warnings, and any assembly and operating instructions. Keep the manual for future reference.

Safety precautions and recommendations detailed here must be fully understood and followed to reduce the risk of injury, fire, explosion, electrical hazard, and/or property damage.

Safety information presented here is generic in nature – some advice may not be applicable to every product. The term "equipment" refers to the product, be it electrical mains powered, battery powered or combustion engine powered.

- **Before Use** - If you are not familiar with the safe operation/handling of the equipment or are in any way unsure of any aspect of suitability or correct use for your application, you should complete training conducted by a person or organization qualified in safe use and operation of this equipment, including fuel/electrical handling and safety.
- Do NOT operate the equipment in flammable or explosive environments, such as in the presence of flammable liquids, gases, or dust. The equipment may create sparks or heat that may ignite flammable substances.
- Keep clear of moving parts.
- Equipment may be a potential source of electric shock or injury if misused.
- Do NOT operate the equipment if it is damaged, malfunctioning or is in an excessively worn state.
- Do NOT allow others to use the equipment unless they have read this manual and are adequately trained.
- Keep packaging away from children - risk of suffocation! Operators must use the equipment correctly. When using the equipment, consider conditions and pay due care to persons and property.

General Work Area Safety

- Work areas should be clean and well lit.
- Do not operate the equipment if bystanders, animals etc are within operating range of the equipment or the general work area.
- If devices are provided for connecting dust extraction / collection facilities, ensure these are connected and used properly. Dust collection can reduce dust-related hazards.

General Personal Safety

- Wear appropriate protective equipment when operating, servicing, or when in the operating area of the equipment to help protect from eye and ear injury, poisoning, burns, cutting and crush injuries. Protective equipment such as safety goggles, respirators, non-slip safety footwear, hard hat, hearing protection etc should be used for appropriate equipment / conditions. Other people nearby should also wear appropriate personal protective equipment. Do not wear loose clothing or jewellery, which can be caught in moving parts. Keep hair and clothing away from the equipment.
- Stay alert and use common sense when operating the equipment. Do not over-reach. Always maintain secure footing and balance.
- Do not use the equipment if tired or under the influence of drugs, alcohol, or medication.
- This equipment is not intended for use by persons with reduced physical, sensory, or mental capabilities.

General Fuel Safety

- Petrol/fuel/gasoline is extremely flammable – keep clear of naked flames or other ignition sources.
- Do not spill fuel. If you spill fuel, wipe it off the equipment immediately – if fuel gets on your clothing, change clothing.
- Do NOT smoke near fuel or when refuelling.
- Always shut off the engine before refuelling.
- Do NOT refuel a hot engine.
- Open the fuel cap carefully to allow any pressure build-up in the tank to release slowly.
- Always refuel in well-ventilated areas.
- Always check for fuel leakage. If fuel leakage is found, do not start, or run the engine until all leaks are fixed.

General Carbon-Monoxide Safety

- Using a combustion engine indoors **CAN KILL IN MINUTES**. Engine exhaust contains carbon-monoxide – a poison you cannot smell or see.
- Use combustion engines **OUTSIDE** only, and far away from windows, doors, and vents.

General Equipment Use and Care

- Keep packaging away from children and pets, potential risk of suffocation or choking.
- The equipment is designed for domestic use only.
- Handle the equipment safely and carefully.
- Before use, inspect the equipment for misalignment or binding of moving parts, loose components, damage, or any other condition that may affect its operation. If damaged, have the equipment repaired by an authorised service centre or technician before use.
- Prevent unintentional starting of the equipment - ensure equipment and power switches are in the OFF position before connecting or moving equipment. Do not carry equipment with hands or fingers touching any controls. Remove any tools or other items that are not a part of the equipment from it before starting or switching on.
- Do not force the equipment. Use the correct equipment for your application. Equipment will perform better and be safer when used within its design and usage parameters.
- Use the equipment and accessories etc. in accordance with these instructions, considering working conditions and the work to be performed. Using the equipment for operations different from those intended could result in hazardous situations.
- Always keep equipment components (engines, hoses, handles, controls, frames, housings, guards etc) and accessories (cutting tools, nozzles, bits etc) properly maintained. Keep the equipment clean and, where applicable, properly lubricated.
- Store the equipment out of reach of children or untrained persons. To avoid burns or fire hazards, let the equipment cool completely before transporting or storing. Never place or store the equipment near flammable materials, combustible gases, or liquids etc.
- The equipment is not weather-proof, and should not be stored in direct sunlight, at high ambient temperatures or locations that are damp or humid.
- Do not clean equipment with solvents, flammable liquids, or harsh abrasives.
- For specific equipment safety use and care, see Equipment Safety.

General Electrical Safety

- Inspect electrical equipment, extension cords, power bars, and electrical fittings for damage or wear before each use. Repair or replace damaged equipment immediately.
- Ensure all power sources conform to equipment voltage requirements and are disconnected before connecting or disconnecting equipment.
- When wiring electrically powered equipment, follow all electrical and safety codes.
- Always use a Residual Current Device (RCD) (GFCI) / Earth Leakage Circuit Breaker / Safety Switch when operating this item (electrical items only).
- High voltage / high current power lines may be present. Use extreme caution to avoid contact or interference with power lines. Electrical shock can be fatal.

General Electrical Safety

- Electrically grounded equipment must have an approved cord and plug and be connected to a grounded electrical outlet.
- Do NOT bypass the ON/OFF switch and operate equipment by connecting and disconnecting the electrical cord.
- Do NOT use equipment that has exposed wiring, damaged switches, covers or guards.
- Do NOT use electrical equipment in wet conditions or in damp locations.
- Do NOT use electrical cords to lift, move or carry equipment.
- Do NOT coil or knot electrical cords, and ensure electrical cords are not trip hazards.

General Service Information

- The equipment must be serviced or repaired at authorised service centres by qualified personnel only.
- Replacement parts must be original equipment manufacturer (OEM) to ensure equipment safety is maintained.
- Do NOT attempt any maintenance or repair work not described in this manual.
- After use, the equipment and components may still be hot – allow the equipment to cool and disconnect spark plugs and/or electrical power sources and/or batteries from it before adjusting, changing accessories, or performing repair or maintenance.
- Do NOT adjust while the equipment is running.
- Perform service-related activities in suitable conditions, such as a workshop.
- Replace worn, damaged or missing warning/safety labels immediately.

Child Choking Hazard!**• WARNING: CHOKING HAZARD - CHILDREN UNDER 8 YRS. CAN CHOKE OR SUFFOCATE ON SMALL PARTS. ADULT SUPERVISION REQUIRED.**

- Carefully inspect anything a baby or child could gain access to, including this product & packaging.
- Dangers include anything of a size that could become stuck or lodged in a baby or child's airway.
- Dangers include things like raw, hard fruit and vegetable pieces, large pieces of meat, bones or sausage skins, popcorn, nuts, hard lollies and corn chips, small magnets and batteries, coins, beads, marbles and small uninflated balloons, broken toys and smaller toys, buttons, keys, and coins. But there are also many more.
- Pay particular attention to accessories, nuts, bolts, screws, washers, caps, covers, loose parts, parts that could be removed or break off, pieces of packaging, staples, tape etc.
- In the case of any doubt, secure the item of concern in an appropriate manner, or completely remove the possibility of access by a baby or child.

Wheelchair Safety

- When attempting to stand you must ensure the footrests have been flipped up prior to attempting to exit the chair.
- Assistance may be required to exit the chair. A qualified carer is recommended when traveling up or down an incline.
- Do not stand directly on the wheels, do not place any weight directly on the wheels.
- When sitting on the seat ensure the brakes are engaged and the user is balanced towards the back of the seat. Do not lean forward.
- The wheelchair is not designed for transiting a patient through long distances or for use on rough, uneven surfaces.

Suitability

Inability to walk or severely pronounced walking disability due to:

- Paralysis
- Loss of limbs
- Limb defect/deformation
- Joint contracture/joint damage (not on both arms)
- Other diseases

Unsuitability

Use of the wheelchair is unsuitable in the case of:

- False sensations
- Severe disequilibrium
- Loss of limbs on both arms
- Joint contracture/joint damage on both arms
- Inability to sit
- Impaired or inadequate vision

Operational Safety





















- Apply the brakes when getting on and off the wheelchair.
- Before getting on or off the wheelchair, do not step on the footplates. They must be folded down beforehand, or the footrests can be moved aside.
- Users requiring additional head stability should only move the wheelchair under the care of an experienced caregiver.

Premium Pediatric Wheelchair

- Carrying the wheelchair with the user in it is unacceptable.
- When moving the wheelchair through doorways and other aisles, make sure that there is enough space on both sides to avoid injury or damage to the wheelchair.
- When using public roads, obey the rules of the road.
- As with other vehicles, you should not drive a wheelchair under the influence of alcohol, medications, or narcotics. This also applies to driving indoors.
- The wheelchair is not intended to transport the user up the stairs, as a stretcher, or by pulling the user on it. The structure of the wheelchair and connecting elements are not adapted to such large loads. If used in this way, it may weaken or damage the entire wheelchair support structure and result in a possible accident.
- The wheelchair is designed for indoor and outdoor use only on hardened and even surfaces.
- Using the product on a ground that does not provide sufficient support to the user, such as icy, wet, oily, or heavily sanded surfaces, is unacceptable—this poses a risk of falling.
- The wheelchair can only be used on surfaces where all four wheels of the wheelchair touch the ground, and contact is sufficient to drive the wheels properly.
- To overcome difficult architectural barriers (steep or long climbs, high steps, etc.), a guardian's help is necessary.
- If the user does not have sufficient experience in using the wheelchair, please ask the guardian for help in overcoming obstacles.
- Do not leave the wheelchair with the user in the immediate vicinity of an open flame.
- Be aware that the surface temperature of the upholstery and frame elements may increase after exposure to the sun or other heat sources. At low temperatures, there is a risk of very high cooling of the metal parts of the wheelchair.
- The wheelchair should be carried by grasping the frame structure of the seat. Never carry the wheelchair by the armrests or footrests.
- The wheelchair cannot be used as a seat in a motor vehicle (bus, train).
- The wheelchair is designed for persons weighing up to 75 kg.
- The wheelchair should always be easy to fold and unfold. In case of problems with the assembly, read the instructions again or contact the manufacturer.
- Be careful when moving on sloping surfaces. The maximum allowable angle of the surface up and down is 5%.
- The brakes are not used to slow down the wheelchair during movement; use them only to prevent unwanted movements of the wheelchair.
- Both brakes should be properly adjusted. Incorrect setting of one brake, resulting in the locking of only one rear wheel on large slopes, may cause difficult-to-control rotation of the wheelchair around the braked wheel.
- The tyre pressure should be checked regularly. If the tyre pressure is lower than recommended, the performance of the brakes can be significantly reduced.
- When lifting objects off the floor, do not reach backwards as you may fall backwards.
- When reaching for objects, avoid excessive leaning out of the wheelchair as it may tip over.
- Do not reach from the front of the floor if it requires sliding out of the wheelchair.
- All activities such as getting on, getting off, changing the user's position, and stopping the wheelchair require the wheel brakes to be engaged.
- Do not leave the user in the wheelchair (even when the brakes are applied) on ramps, slopes, or uneven surfaces.
- After assembly of the wheels, check whether the locking of the connecting axle of the wheel has engaged. During assembly, you should hear a characteristic click of the latch. In addition, strongly pull the wheel without pressing to ensure it is secure.
- Using the parking brake while driving may cause an accident.
- Always check after each adjustment that the brake works correctly.
- In the event of contamination of the tyres by oils, grease, water, or other substances, braking efficiency may be reduced.
- After prolonged use of the wheelchair, the braking force may weaken due to wear between the brake pin and the tyre of the wheel. Test the effectiveness of the brake adjustment on a wheelchair with a user.
- To maintain the wheelchair's operational parameters, it is extremely important to keep all its parts in a safe place and ensure proper storage.
- Storing the wheelchair in rooms with very high humidity (e.g., bathrooms, laundries) or in the open air, combined with a lack of hygiene, will cause a rapid loss of performance and aesthetic quality of the wheelchair.
- It is unacceptable to use the wheelchair if any visible or noticeable damage has been found. Any wear of the wheelchair's parts indicates a potential risk of accident (e.g., worn wheels, upholstery cracks, etc.). All repairs should be entrusted to qualified personnel, and during the warranty period, the manufacturer's service centre.
- In the case of wheelchairs equipped with inflatable wheels, do not exceed the maximum permissible tyre inflation pressure as stated on the tyre.
- During transport (e.g., in a luggage carrier), do not load the wheelchair with anything, as it may affect its structural integrity.
- Particular care should be taken when folding the wheelchair. Lack of attention may result in fingers being trapped during the following actions:
 - Putting on and removing the footrest.
 - Unfolding the wheelchair within the moving parts in the brakes area.
- Keep the wheelchair out of reach of children under five years of age.

Safety Symbols

The product may have safety warning labels attached to it, explained below. Understand the symbols on your product and their meanings. If any stickers become unreadable, unattached etc., replace them.

 <p>Flammable Material Hazard Flammable liquids, gases, or substances etc may present. Avoid ignition sources and open flames. Danger of fire.</p>	 <p>Read User Manual Read and fully understand product safety warnings, operation, procedures etc before using the product.</p>	 <p>Use Hand Protection Wear appropriate hand protection and take due care as the product or use of the product may present hand hazards.</p>	 <p>Carbon-Monoxide Hazard Do not use the product in confined areas or without adequate ventilation. Carbon-monoxide poisoning can be fatal.</p>
 <p>Electrocution / Electrical Shock Hazard High voltage or high current electricity may be present or required by the product. Take due care when handling electrical products, cables, plugs and leads. Electrical shock can be fatal.</p>	 <p>Toxic Fumes / Dust Hazard Using the product or by-products from use may produce fumes, smoke or particles that could be harmful if inhaled. Wear appropriate breathing protection and have adequate ventilation.</p>	 <p>Explosive Material Hazard Combustible liquids, gases, or substances etc may be present. Avoid ignition sources and open flames. Danger of explosion.</p>	 <p>Cutting / Amputation Hazard The product may have blades, edges or mechanical devices that can cause severe cut injury to fingers, limbs etc. Take due care when handling and using the product.</p>
 <p>Crush Hazard The product may have blades, edges or mechanical devices that can cause severe crush injury to fingers, limbs etc. Take due care when handling and using the product.</p>	 <p>Single Operator Only The product must be operated by a single person only. More than one person operating the product may introduce additional hazards.</p>	 <p>Use Face Protection Wear appropriate full-face protection and take due care as the product or use of the product may present face and eye hazards.</p>	 <p>Use Foot Protection Wear appropriate foot protection and take due care as the product or use of the product may present foot hazards.</p>
 <p>Use Eye / Ear / Head Protection Wear appropriate eye and / or ear and / or head protection and take due care as the product or use of the product may present eye, hearing, and head hazards.</p>	 <p>Running Hazard Do not run on or near the product as doing so may present a fall hazard.</p>	 <p>Diving Hazard Do not dive into the product as doing so may present a neck / head injury hazard.</p>	 <p>Adult Supervision Required Always supervise children and other users of a product to prevent drowning or injury.</p>
 <p>Skin Penetration / Puncture Hazard The product may produce pressure, emit liquids or objects that can cause severe injury to fingers, limbs, blood etc. Take due care when handling and using the product.</p>	 <p>Hot Surface Hazard Be aware that the product may produce high temperatures and hot surfaces that can cause burn injuries.</p>	 <p>Flying Debris Hazard Be aware that the product or use of the product may present hazards produced by flying debris. Wear appropriate clothing and protective devices.</p>	 <p>Moving Parts Hazard Be aware that the product contains or uses mechanical devices that move or rotate. Always wait for moving parts to stop fully before handling the product, adjusting, maintenance etc.</p>
















 <p>Carbon-Monoxide Hazard Do not use the product in confined areas or without adequate ventilation. Carbon-monoxide poisoning can be fatal.</p>	 <p>Pull Hazard Be aware that the product contains or uses mechanical devices that can pull in objects and can cause severe injury to fingers, limbs etc. Take due care when handling and using the product.</p>	 <p>Slope / Fall Injury Hazard Be aware that using the product on sloping surfaces or in slippery conditions may present additional dangers from falls and contact with blades, moving parts, hot surfaces etc.</p>	 <p>"Slam Dunk" Warning Do NOT attempt "slam dunk" manoeuvres as this may result in severe injury due to falling, product breakage or collapse etc.</p>
 <p>Electrocution / Electrical Shock Hazard - Outdoor High voltage or high current electricity may be present or required by the product. Do NOT use in rain, damp, or wet conditions. Electrical shock can be fatal.</p>	 <p>Electrocution / Electrical Shock Hazard - Disconnect High voltage or high current electricity may be present or required by the product. Always disconnect the product from the electrical supply before handling the product, adjusting, maintenance etc.</p>	 <p>Power Line Electrocution Hazard High voltage / high current power lines may be present. Use extreme caution to avoid contact or interference with power lines. Electrical shock can be fatal.</p>	 <p>"Kick-Back" Hazard High level of "kick-back" hazard that can cause the machine to suddenly rotate towards operator. Kick-back injury can be fatal.</p>
 <p>Winch Operator Position Hazard Do NOT stand between winch and load. Do NOT use winch to move people.</p>	 <p>Winch Lift Hazard Do NOT LIFT load vertically. Use machine to PULL only.</p>	 <p>Cable Hazard Ensure that load bearing cable is not kinked or knotted.</p>	 <p>Winch Cable Hazard Ensure that there is a minimum number of cable coils on winching mechanism.</p>
 <p>Winch Hook Hazard Carry hook to load – do NOT throw or run.</p>	 <p>Flash / Blinding Hazard Wear appropriate eye protection for welding. Direct exposure to weld arcs may cause permanent eye injury.</p>	 <p>Laser Hazard Laser may be in use – do NOT look directly at laser or allow others to.</p>	 <p>Choking Hazard Children under 8 yrs. Can choke or suffocate on small parts. Adult supervision required.</p>

Table of Contents

Safety	2
Safety Symbols	5
Parts Identification	8
Assembly and Adjustments	9
Unfolding the Wheelchair	9
Folding the Wheelchair	9
Installation Instructions	9
Adjusting the Height of the Armrest	9
Folding the Footplates Up	10
Swivelling Away and Detaching the Leg Rests	10
Propelling and Slowing Down the Wheelchair with the Hand Rims (If the Rear Wheel Has Hand Rims)	10
Hand Brake (Locking Brake or Tyre Brake)	10
Braking by a Caregiver (Attendant Brake)	10
Removing the Rear Wheel	11
Maintenance	12
General	12
Tyre Pressure	13
Specifications	14

Parts Identification

The wheelchair is packed in a cardboard box. Please keep the packaging if possible; it can be used to store the wheelchair if necessary.



NOTE: Colours and designs may vary

Assembly and Adjustments

IMPORTANT! When attempting to stand you must ensure the footrests have been flipped up prior to attempting to exit the chair. Assistance may be required to exit the chair. A qualified carer is recommended when traveling up or down an incline. Do not stand directly on the wheels, do not place any weight directly on the wheels. When sitting on the seat ensure the brakes are engaged and the user is balanced towards the back of the seat. Do not lean forward. The wheelchair is not designed for transiting a patient through long distances or for use on rough, uneven surfaces.

Unfolding the Wheelchair

1. Stand in front of the wheelchair frame.
2. Push down and outward on the sides of the seat frame, keeping your fingers pointing inwards to prevent them from getting pinched.
3. Continue pushing downwards until the wheelchair frame is fully opened.
4. Place the cushion on the seat and place the backrest cushion on the backrest tube.
5. Fix the headrest, rib external fixator, leg splitter, calf strap, and footplates.
6. If needed, install the table on the side of the seat.

Folding the Wheelchair

1. Remove cushions, chest safety belt, and calf strap, then fold the footplates.
2. Stand beside the wheelchair frame, grasp the tubes of the seat, and pull upwards.

WARNING: Never attempt to lift the wheelchair by any detachable parts (e.g., armrests). This may result in injury.

Installation Instructions

1. Open the wheelchair frame, install cushions, chest safety belt, and leg rests.
2. Unfold the footplates and install the calf strap.

Adjusting the Height of the Armrest

The armrest can be removed, and the height can be adjusted. There are three height options to choose from.

1. Hold the handle of the height adjustment key, lift it up, then adjust the height of the armrest.
2. Finally, hold the handle of the height adjustment key, and pull it downward.




Folding the Footplates Up

The footplates can be removed using a wrench to fix or remove them. Fold the footplates up to make entering and exiting easier. The footplates must be folded up to save space when folding the wheelchair.

Swivelling Away and Detaching the Leg Rests

1. To remove the leg rests, first fold the footplates up.
2. Then, pull the locking lever and swivel the leg support outwards.
3. The leg support can now be detached upwards.
4. Carry out the procedure in the opposite sequence during assembly. When the leg support is swivelled forwards, you should hear it engage clearly.



Propelling and Slowing Down the Wheelchair with the Hand Rims (If the Rear Wheel Has Hand Rims)

- The wheelchair can be manoeuvred forwards, backwards, and left to right by using the hand rims on the rear wheels.
- Hold the hand rim with your hand and move the wheel forwards or backwards.
- Your thumb should rest on top while you grab the hand rim by making a fist.
- To slow down the wheelchair, control the hand rim by letting it slide through your hand.

Hand Brake (Locking Brake or Tyre Brake)

- The hand brake enables the user to park the wheelchair securely.
- Actuate the hand brake (tyre brake) when stationary: To do this, push the lever forwards until it stops independently. The wheelchair is now securely parked; you can transfer or simply leave it parked.
- The hand brake (locking brake) is released by pulling the lever backwards.
- While moving, always slow down the wheelchair by using the hand rims. The hand brakes are only used for parking and cannot control the wheelchair at speed.

CAUTION! When propelling your wheelchair, make sure that you do not touch the tyre casing with your thumb. Risk of trapping between the tyre and side panel/armrest!

Braking by a Caregiver (Attendant Brake)

If the wheelchair has the attendant brake, the caregiver can help the user stop or slow down the wheelchair by pulling both attendant brakes upwards.

Note: Always use the left and right brakes simultaneously. Using only one side will make the chair very hard to manoeuvre.

Removing the Rear Wheel

Push the quick-release button on the axle, and the rear wheel can be removed.



Maintenance

IMPORTANT! Ensure the wheel axles are checked and tightened regularly.

General

Inspections should be carried out regularly and should be preceded by a thorough cleaning of the wheelchair (to reveal any hidden damage) during the technical inspection.

Frequency	Subject of Control	Comments
Before First Assembly	Checking the technical condition of all elements contained in the packaging	Check that no components have visible mechanical damage.
Before Each Ride	Wheel and brake control	The brakes may wear out; check the condition of the brakes before each use. Check the wear of the wheels and clean them of dust and sand.
Every Two Weeks	Fastening screws and nuts for correct mounting	Ensure that all nuts and bolts are tight.
Every Two Months	Cleanliness, overall condition	In the event of mechanical damage, contact the manufacturer immediately.
Once a Month	Wheel axles	The wheel axles should be maintained by applying a few drops of lubricant between the axle and the hub. The wheelchair must be lifted so that the wheel under test is above the ground. The wheel should turn freely but without slots on the axle, and the front wheels should also rotate freely on the axis of direction.
Once a Month	Condition of quick-connect, threaded, and soldered connections	Pay attention to cracks in the frame material at soldering points and at the holes blocking the bayonet connections.
Once a Month	Tyre condition	Ensure that the tyre tread is not worn and that the tyres have no visible damage.

Tyre Pressure



The tyre pressure should be checked every 30 days, If the tyre is going flat, unscrew the cap on the Schrader valve and inflate the tyre to the required psi. The required psi is typically printed on the sidewall of the wheelchair's tyres.

**Not applicable to solid rubber tyres.*

Specifications

Chassis Material	Powder-coated Aluminium
Frame & Chassis Finish	Powder Coated
Max User Load	75kg
Seat Width	36cm
Seat Height	53cm (from floor)
Brake Type	Park brakes
Rear Wheels	22" Pneumatic tyres
Front Wheels	6" Never-flat front tyres
Supports	Single crossbar
Seatbelt	Yes

Note: Requires minimal assembly



Some experts believe that the incorrect or prolonged use of almost any product may cause serious injury or death. To help reduce your risk of serious injury or death, refer to the information below. For more information, see www.datastreamserver.com/safety

- Consult all documentation, packaging, and product labelling before use. Note that some products feature documentation available online. It is recommended to print and retain the documentation.
- Before each use, check the product for loose/broken/damaged/missing parts, wear, or leaks (if applicable). Never use a product with loose/broken/damaged/missing parts, wear, or leaks.
- Products must be inspected and serviced (if applicable) by a qualified technician every 6 months. This is based on average residential use by persons of average size and strength, and on a property of average metropolitan size. Use beyond these recommendations may require more frequent inspections/servicing.
- Ensure that all users of the product have completed a suitable industry recognised training course before being allowed access to the product.
- The product has been supplied by a general merchandise retailer that may not be familiar with your specific application or description of application. Be sure to attain third-party approval from a qualified specialist for your application before use, regardless of any assurances from the retailer or its representatives.
- This product is not intended for use where fail-safe operation is required. As with any product (for example, automobile, computer, toaster), there is the possibility of technical issues that may require the repair or replacement of parts, or the product itself. If the possibility of such failure and the associated time it may take to rectify could in any way inconvenience the user, business, or employee, or financially affect the user, business, or employee, then the product is not suitable for your requirements. This product is not intended for use where incorrect operation or a failure of any kind, including but not limited to, a condition requiring product return, replacement, parts replacement, or service by a technician may cause financial loss, loss of employee time or an inconvenience requiring compensation.
- If this product has been purchased in error when considering the information presented here, contact the retailer directly for details of their returns policy, if required.

