## **PHILIPS**

# Star Wars shaver

## Wet and dry electric shaver

ComfortCut Blade System

5-direction Flex Heads

Turbo+ mode

Cleansing brush & P. trimmer



#### SW5720/46







### Close shave, smooth skin

### Wield the force!

Seek empowerment through shaving innovation. Our ComfortCut blades with rounded profile heads glide smoothly over your skin for a protective shave.

#### A comfortable and close end result

- $\boldsymbol{\cdot}$  Rounded edges move smoothly over skin for a protective shave
- · 27 rotating blades capture and cut hair from all angles
- ${\boldsymbol{\cdot}}{}$  Heads flex in 5 directions for a fast close shave
- Shave through dense parts of beard with 20% more power
- Get a comfortable dry or refreshing wet shave with Aquatec

#### Easy to use

- 70 minutes of cordless shaving
- 1-hour charging time
- $\boldsymbol{\cdot}\,$  Shaver can be rinsed clean under the tap
- · 2-year guarantee, worldwide voltage and replaceable blades

#### Get the most out of your shaver

- $\boldsymbol{\cdot}\;$  Click-on brush to cleanse more thoroughly than by hand
- Click-on trimmer for perfect mustache and sideburn trimming



### Highlights

#### **ComfortCut Blade System**



Get a comfortable shave, wet or dry. Our ComfortCut blades with rounded profile heads glide smoothly across your skin while protecting it from nicks and cuts.

#### 27 self-sharpening blades



27 self-sharpening blades. 56000 cuts per minute. Zero hairs left standing – no matter what direction they're growing in.

#### 5-direction Flex Heads




5-direction Flex Heads with 5 independent movements ensure close skin contact for a fast and close shave even on neck and jaw line

#### **Turbo+ mode**



Shave even faster through dense parts of your beard with the 20% extra power boost you get by activating Turbo+ mode.

#### **Aquatec Wet & Dry**



Choose how you prefer to shave. With the Aquatec Wet & Dry seal, you can opt for a quick yet comfortable dry shave. Or you can shave wet — with gel or foam — even under the shower.

#### 70 minutes of cordless shaving



You'll have 70+ minutes of running time – that's about 25 shaves – on a one-hour charge. This shaver only operates in cordless mode.

#### 1-hour charging time



Shave longer with every charge, thanks to our powerful energy-efficient, long-lasting lithium-ion battery. Do a quick charge for 5 minutes and you'll have enough power for one shave

#### SmartClick cleansing brush



Click on our cleansing brush to clean your face more thoroughly than you can by hand. Use it with the cleanser of your choice. It removes dead skin cells and increases microcirculation, leaving your skin feeling fresh and radiant.

#### **SmartClick precision trimmer**



Click on our skin-friendly precision trimmer to finish your look. It's ideal for maintaining your mustache and trimming your sideburns.









### Specifications

#### **Shaving Performance**

Shaving system: ComfortCut Blade System, 27

self-sharpening blades

Contour-following: 5-direction Flex Heads Skin comfort: Skin Protection System,

AquaTec Wet & Dry

Turbo+ mode: Shave with 20% extra power

#### **Accessories**

SmartClick: Deep facial cleansing brush, Precision trimmer Pouch: Travel pouch

#### Power

Run time: 70 min / 25 shaves

Charging: 1 hour full charge, Quick charge 5

min for 1 shave

Battery Type: Lithium-ion Automatic voltage: 100-240 V Max power consumption: 9 W Stand-by power: 0.1 W

#### Ease of use

Wet & Dry: Wet and dry use Cleaning: Fully washable

Operation: Cordless use, Unplug before use

#### Design

Handle: Ergonomic grip & handling

Color: Star Wars Kylo Ren

#### **Service**

2-year guarantee

Replacement head: Replace every 2 yrs with

SH50

© 2017 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2017-08-23 Version: 6.2.1



<sup>\*</sup> SmartClick cleansing brush - compared to manual cleansing in same condition, according to the results from the external clinical test on 30 candidates and in 3 hours after cleansing